

● cDNA -20

RT-PCR Mix for qPCR

RT-PCR Enzyme Mix		DNase	DNA
RNA	5X RT-PCR Buffer	Control Mix	
RNA	DNA		
PCR			

5X RT-PCR Buffer	4 µl
RT-PCR Enzyme Mix	1 µl
RNA	1pg-1µg
RNase-free ddH ₂ O	20 µl

RT-PCR Enzyme Mix	100 µl
5X RT-PCR Buffer	400 µl
Control Mix	20 µl
RNase-free ddH ₂ O	1 ml × 2

-20

cDNA

- DNase
- RT-PCR Enzyme Mix 5X RT-PCR Buffer Control Mix
- RNase RNA RNA -70

qPCR RNA

DNA

5X RT-PCR Buffer	4 µl
Control Mix	1 µl
RNA	1pg-1µg
RNase-free ddH ₂ O	20 µl

PCR

50	15 min
85	5 sec

cDNA qPCR cDNA

