



3-5 μ l
500bp Ladder
8 cm 1 %
1 $\times$ TAE 7 V/cm 40 min
4 2

500bp Ladder
M1251/60 M1252/60 3
500bp Ladder 10 DNA DNA DNA Marker/ DNA
Ladder 500 bp 5 kb DNA
92 ng/ μ l

500 1,000 1,500 2,500 3,000 3,500 4,000 4,500 5,000

- 1.
- 2. DNA Marker
- 3. 0.75-1mm 2-3 μ l Marker 5mm 0.75-1.5mm 3mm
- 4. 4-6 μ l Marker
- 5. PAGE 0.5-1 μ l Marker 1 loading buffer

Marker Lambda DNA Marker DNA
DNA DNA DNA PCR
Marker DNA Marker DNA DNA Marker
1 SensiAnsys DNA
DNA Marker 100 bp Marker
/