




3-5 μ l/
10bp Ladder

8 cm 5 %
0.5×TBE 5 V/cm 1 h
20 cm 10 %
1×TBE 8 V/cm 3 h

3

-20

100 ng/5μl
40 ng/5μl

10bp Ladder

M1011/60	M1012/60	3
----------	----------	---

10bp Ladder	18	DNA
Ladder	80 bp	300 bp

168 ng/ μ l[illegible]

1. DNA Marker

2. DNA Marker

3. 0.75-1mm 2-3 μ l Marker 5mm 0.75-1.5mm 3mm

4. Marker

Marker

Lambda DNA Marker

DNA

DNA DNA DNA PCR

Marker
DNA Marker DNA
 DNA DNA Marker

SensiAnsys DNA

1

DNA Marker
100 bp

1

Marker