


00



5 μl/
LD 500bp Ladder

8 cm 1 %
1×TAE 7 V/cm 40 min

4 2

12.5 ng/5μl
5 ng/5μl

Ladder

LM1252/60 3

Ladder

10
500 bp

5 kb

DNA

DNA

DNA

DNA Marker/ DNA

2000

11.5 ng/μl

500 1,000 1,500 2,000 2,500 3,000 3,500 4,000 4,500 5,000

1. Marker
2. DNA
3. Marker 3mm
0.75-1mm 2-3 μl 5mm 0.75-1.5mm
4. Marker

5. PAGE 0.5-1 μl Marker 1 loading buffer

1. Marker Lambda DNA Marker DNA

2. DNA DNA DNA PCR

3. DNA

4. Marker DNA Marker DNA DNA Marker

1 SensiAnsys DNA

/ DNA Marker 100 bp Marker